BPBIO320S

User's Manual

BPBIO320S User Manual for Measurement Guide and Setup

Thank you for purchasing the BPBIO320S. This user's manual describes all the features of the BPBIO320S. Please read before use and keep it in a safe place. By following the manual instructions, you will be able to use the BPBIO320S more safely and effectively.

Please note the important information below before reading this manual.

Warning

Failure to observe these precautions can result in personal injury or equipment damage.

Caution

Failure to comply with safety precautions can damage the equipment.

Note

Referring to notes can help improve equipment use.

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InBody

For problems with the BPBIO320S or any clinical questions, contact:

Representative & Sponsor Information

EC REP InBody Europe B.V.

Gyroscoopweg 122, 1042 AZ, Amsterdam, The Netherlands

TEL: +31-20-238-6080 FAX: +31-6-5734-1858 Website: nl.inbody.com E-mail: info.eu@inbody.com

Australian Sponsor: [AUSTRALIA]

Emergo AUSTRALIA. Level 20, Tower II, Darling Park, 201 Sussex Street, Sydney, NSW 2000, AUSTRALIA

InBody Co., Ltd.

15, Heugam-gil, Ipjang-myeon, Seobuk-gu, Cheonan-si, Chungcheongnam-do 31025 KOREA

TEL: +82-41-581-3003 FAX: +82-41-581-3103 Website: inbody.com E-mail: info@inbody.com

Intended Use

The BPBIO320S is a digital monitor intended for use in measuring blood pressure and pulse rate in user population with upper arm circumference ranging from 17cm to 42cm (7-inch to 17-inch). The systolic blood pressure and diastolic blood pressure are measured by non-invasive blood pressure (“NIBP”) measuring method. The BPBIO320S may provide useful clinical information about the current health status of not only the users who are diagnosed with hypertension but also those who are not diagnosed with hypertension. Warnings and cautions described in the user’s manual should be observed at all times.

*This device is not designed for use in home environment.

Customer Service Information

InBody USA [USA]

13850 Cerritos Corporate Dr. Unit C Cerritos, CA 90703 USA

TEL: +1-323-932-6503 FAX: +1-323-952-5009 Website: inbodyusa.com E-mail: info.us@inbody.com

InBody Japan [JAPAN]

Tani Bldg., 1-28-6, Kameido, Koto-ku, Tokyo 136-0071 Japan

TEL: +81-3-5875-5780 FAX: +81-3-5875-5781 Website: www.inbody.co.jp E-mail: inbody@inbody.co.jp

InBody China [CHINA]

903/904, XingDiPlaza, No.1698 YiShanRoad, Shanghai 201103 China

TEL: +86-21-64439705 FAX: +86-21-64439706 Website: inbodychina.com E-mail: info@inbodychina.com

InBody Asia [ASIA]

Unit 3A-11, Oval Damansara, 685 Jalan Damansara Kuala Lumpur, WP KL 60000 Malaysia

TEL : +60-3-7732-0790 FAX: +60-3-7733-0790 Website: inbodyasia.com E-mail: info@inbodyasia.com

InBody India [INDIA]

InBody India Pvt. Ltd.

57/57 A, 1st Floor, Raj Industrial Complex, Military Road, Marol, Andheri (East). Mumbai- 400059, Maharashtra, India

TEL : +91-22-6223-1911 Website: inbody.in E-mail: india@inbody.com

BPBIO320S

User's Manual Contents

I. BPBIO320S Installation

A. Product Components	6
B. Operating Environment	6
C. Installation Instructions	7
D. Loading and Reloading Printer Paper	11
E. Replacing Cuff Cover	13
F. Setup	15
G. Connecting Multi-Display Device	21
H. Connecting to InBody	22
I. Precautions for Maintenance	23

II. Blood Pressure Test

A. Precautionary Steps	24
B. Test Instructions	24
C. Test Posture	28

III. Transportation and Storage

A. Cautions during Transportation	29
B. Repacking Instructions	30
C. Transportation and Storage Environment	30

IV. FAQ

A. Error Code	31
B. Regarding the Equipment	31
C. Regarding the Test	34

V. Others

A. Exterior and Functions	36
B. Safety Information	37
C. Product Classification	38
D. Specifications	39
E. EMC Declaration	40

I. BPBIO320S Installation

A. Product Components

The BPBIO320S consists of the following components. Make sure the following components are present.


* Please check each component of the system for damage prior to installation.

Basic Components

- ❶ BPBIO320S
- ❷ Power cable (AC 250V 10A 1.8m) 1 EA
- ❸ Information Panel
- ❹ Two rolls of printing paper
- ❺ Two cuff covers (One installed)
- ❻ User's Manual

Optional components

- ❼ Multi-Display Device
- ❽ Multi-Display Device Cable
- ❾ Desk for BPBIO320S,
and two screws
- ❿ Chair for BPBIO320S


B. Operating Environment

Please make sure the environment is adequate for the BPBIO320S installation. This equipment is designed for indoor use. If installing the equipment outdoors, the following requirements must be fulfilled:

Temperature range	5 to 40°C (41 to 104°F)
Relative humidity	15 to 90% RH
Atmospheric pressure range	70 to 106 kPa

C. Installation Instructions

1. Open the BPBIO320S box, and remove the air cushion (small size). Unfold the air cushion (large size) and take out the BPBIO320S from the box, and then remove the protective box. Take off the protective vinyl out of the BPBIO320S, and remove the packing material.


! Caution

- Do not carry the equipment by holding the cuff.
- Keep packing material to repack the BPBIO320S at a later time. Dispose of trash according to local laws and regulations.
- Refer to the figure below when transporting the BPBIO320S.

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
- Place the equipment on a solid table that can withstand the weight of the BPBIO320S (approx. 9kg (20 lb)). If desk for the BPBIO320S is available, align the grooves on the bottom of the BPBIO320S with the grooves on the desk and secure them with two fixing screws.

* Desk for the BPBIO320S sold separately.


! Caution

- Install the BPBIO320S on a flat, vibration-free floor. If the equipment is installed where the floor is not flat, it may topple during a test or the test results may be inaccurate.

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
- Mount the Information Panel onto the BPBIO320S.

- Connect the power cable provided with the BPBIO320S. Plug the power cable into a grounded 3-prong outlet.

* The anti-theft groove can be used to secure the equipment to the table.

Caution

- If you are not using a power plug type B or F, you must use a 3-prong plug type with a ground pin.

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Warning

- Do not install the equipment where the power cable is difficult to disconnect.
- Do not plug or disconnect the power cable with wet hands. There is a risk of electric shock.
- Always use an outlet supplied to the correct power supply (AC 100-240 V). An incorrect power supply may cause fire or malfunctions.
- When using a multiple tab, use it having enough power capacity or use an extension cable.
- Do not disassemble the equipment arbitrarily. This may cause electric shock, injury, product malfunction, and/or inaccurate results. Unauthorized disassembly will void the manufacturer's warranty.
- Do not let the BPBIO320S touch with other electronic devices when on. This may cause electric shock.

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Caution


- Plugging the device into an ungrounded outlet may cause malfunction or electric shock.
- Always handle the equipment with care. Exterior impact may cause trouble when operating the equipment.

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
5. Use the anti-theft groove to prevent theft.


* Theft protection cable is sold separately.

6. Turn the BPBIO320S on by pressing the power switch located on the back of the device.

7. When the BPBIO320S is turned on, the following screen is displayed and a calibration is automatically performed.


Caution

- Do not operate or touch the device while it is calibrating. It may cause errors to future tests.

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
8. Adjust the volume by turning the volume knob on the rear of the BPBIO320S.


D. Loading and Reloading Printer Paper

1. Loading Printer Paper

- 1) Cut off the end of the printer paper with scissors. If the end of the paper is ripped, it may cause a paper jam or damage the printer head.


- 2) Open the cover by pressing 'PUSH' while the BPBIO320S is turned on.


- 3) If the device is being used for the first time, remove the protective paper located inside.

* The protective paper (thermal printing paper) is in the printer groove. Pull upward to remove it.

- 4) Fold one corner of the provided printer paper and insert it in the direction as shown.


Caution

- The printer will not work if the paper is inserted incorrectly.

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5) Insert the end of the paper through the feeder. The printer will automatically pull and cut the paper.


6) Remove the cut paper, and close the printer paper cover.

2. Reloading Printer Paper

1) Open the cover by pressing 'PUSH' mark on the bottom center of it while the BPBIO320S is turned on.

2) Press  button to remove the remaining paper.

3) Reload a new roll of printer paper by following the instructions for 'Loading Printer Paper'.

Caution

- Remove all trash from previous printer paper rolls.

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Note

- Two rolls of printer paper are provided with the BPBIO320S.
- If additional printer paper roll is needed, it must be purchased separately as it is a consumable item.

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

E. Replacing Cuff Cover

The cuff cover of the BPBIO320S should be replaced periodically for sanitary purposes. Follow the guide below to replace the cuff cover.


1. Remove the front cover of the BPBIO320S by pushing downward. Loosen the screw with a screwdriver.


2. Loosen the two screws on the rear of the BPBIO320S with a screwdriver.

3. Remove the rings on the front and rear of the BPBIO320S by pulling the arrow mark and pull the cuff cover to remove it.

4. Place the seam of the new cuff cover down as shown, and align it with the front and rear ring.


! Caution

- Blood pressure results may be inaccurate if the seam is not face down.

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

5. Secure the new cuff cover onto the grooves on the front and rear of the BPBIO320S.

6. Attach the previously removed front and rear rings on the BPBIO320S and secure them with screws.

7. Close the front cover of the BPBIO320S.

Note

- Two cuff covers are provided with the BPBIO320S. (One is already mounted in the equipment.)
- If additional cuff cloth is necessary, the user must purchase it separately as a consumable item.

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F. Setup

Setup consists of sound settings, motion sensor settings, elbow sensor settings, Automatic printing options, time stamps and Result sheet item(Outputting caution comment, PULSE graph, Additional result values).

Note


- Contents of results sheet output are only applicable to BPBIO320S.

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
1. Settings

- 1) When  button and  button are pressed and held at the same time for 2 seconds while on the home screen, the first screen of the Setup, 'Sound Setting' is displayed.


- 2) Use UP (▲) and DOWN (▼) buttons to change the settings.

* UP (▲) and DOWN (▼) are located beneath the cover of the device.

- 3) Press  button to move to the next screen. Set values are automatically saved.

- 4) Press and hold  button for 2 seconds to return to the home screen.

* The measurement value display time setting is applied when the blood pressure monitor restarts.

※ The following is an example showing how to change a 'manual output' Automatic printing options.

① While on the home screen, press and hold  button and  button at the same time for 2 seconds.

② Press  button to move to *F04Prt* screen.

③ Press  button to select *OFF*.

④ Press and hold  button for 2 seconds to save changes and return to the home screen.

※ The following is an example showing how to change the screen timeout duration.

① While on the home screen, press and hold  button and  button at the same time for 2 seconds.

② Press  button several times to move to *F10Lcd* screen.

③ Press  button to select *10*.

④ Press and hold  button for 2 seconds to save changes and return to the home screen.

2. Detailed description of Setup

When in the Setup menu, the following information is displayed in the Systolic, Diastolic and Pulse sections of the screen.


Function	Display			Description
	Systolic	Diastolic	Pulse	
Sound	F01	Snd	U01	Reading test instruction and measurement results
			U02	Reading test instruction only
			U03	Reading measurement results only
			bEF	Beep guide
			OFF	Mute
Setting up the motion sensor	F02	HI-	On1	Motion sensor is operating with “ding-dong” sound
			On2	Motion sensor is operating in silence
			OFF	Motion sensor is disabled
Setting up the Elbow Sensor	F03	Arn	On	Elbow sensor is operating
			OFF	Elbow sensor is disabled
Outputting results sheet	F04	Prt	On	Automatic printing of result sheet
			OFF	Manual printing of result sheet
Hour setting	F05	HOU		setting hour (0-23)
Minute setting	F06	Min		setting minute (0-59)
Year setting	F07	YER		setting year (2000-2063)
Month setting	F08	Mon		setting month (1-12)
Day setting	F09	Day		setting day (1-31)
Screen timeout duration	F10	LCD	5	Display the test result value for 5 seconds
			10	Display the test result value for 10 seconds
			20	Display the test result value for 20 seconds
			179	Display the test result value for 2 minutes

Result sheet item	Caution comment	F 11	6A	SIC	Output of caution comment
			SCH	00L	Not used
	PULSE graph	F 12	9rP	On	Output of pulse graph
				OFF	Not used
	Additional result values	F 13	0Pt	On	Mean arterial pressure/PULSE pressure /Pressure rate product output
				OFF	Not used
	QR code*	F 14	9r	On	Output of QR code*
				OFF	Not used
Measurement mode setting	3 consecutive measurements	F 15	Avg	On	3 consecutive measurements mode
				OFF	One time measurement mode


* “QR Code” is registered trademark of DENSO WAVE INCORPORATED.

3. Blood Pressure Evaluation

- 1) While on the home screen, press and hold  button and UP (▲) button at the same time for more than 2 seconds.


- 2) When 111 is displayed in the systolic blood pressure column, diastolic blood pressure column, and pulse column as shown below, the blood pressure evaluation has been disabled.


- 3) When 222 is displayed in the systolic blood pressure column, diastolic blood pressure column, and pulse column as shown below, the blood pressure evaluation has been turned on.

* Blood pressure evaluation comment type: Normal, Prehypertension, Stage 1 Hypertension, Stage 2 Hypertension

2017/12/06 19:12
SYSTOLIC 110 mmHg
DIASTOLIC 75 mmHg
PULSE 96 bpm
* Test result - Normal Please consult with a physician for diagnosis and treatment.


PULSE Graph(F12): OFF
Option(F13): OFF

2017/12/06 19:12
SYSTOLIC 110 mmHg
DIASTOLIC 75 mmHg
PULSE 96 bpm
PULSE GRAPH 
* Test result - Normal Please consult with a physician for diagnosis and treatment.

PULSE Graph(F12): ON
Option(F13): OFF

2017/12/06 19:12
SYSTOLIC 110 mmHg
DIASTOLIC 75 mmHg
M.A.P 87 mmHg
P.P 35 mmHg
PULSE 96 bpm
R.P.P 10560
* Test result - Normal Please consult with a physician for diagnosis and treatment.

PULSE Graph(F12): OFF
Option(F13): ON



2017/12/06 19:12
SYSTOLIC 110 mmHg
DIASTOLIC 75 mmHg
M.A.P 87 mmHg
P.P 35 mmHg
PULSE 96 bpm
R.P.P 10560
PULSE GRAPH 
* Test result - Normal Please consult with a physician for diagnosis and treatment.

PULSE Graph(F12): ON
Option(F13): ON

- 4) In order to change the blood pressure evaluation settings, press and hold  button and UP (▲) button at the same time for more than 2 seconds.

4. Displaying test counts

- 1) While on the home screen, hold the DOWN (▼) button for more than 2 seconds. Test counts for the day will be displayed in the diastolic blood pressure column. Test counts for the day and total test counts are printed.
- 2) When DOWN (▼) button is pressed again, the total test counts are displayed in the diastolic blood pressure column.

2017/12/06 09:12
* Today Count : 0002
* Total Count : 00075

Test count printout


G. Connecting Multi-Display Device

Check the results of a blood pressure test by connecting the multi-display device to the BPBIO320S.


* Multi-display device is optional, it can be purchased separately.

* The multi-display device is packaged in a separate air-cell in the BPBIO320S packaging box.


1. Prepare a multi-display device and a multi-display device cable.


2. Connect the multi-display device cable to the SUB DISPLAY connection terminal on the back of the device while the BPBIO320S is turned off. Then connect the other end of the display device cable to connection terminal on the left of the display device.

3. Insert the connected multi-display device into the groove on the top of the BPBIO320S. Then, insert the Information Panel into the groove on top of the display device.

*The BPBIO320S can be used as as stand alone device as shown below.

H. Connecting to InBody

If the BPBIO320S connects to InBody products, the blood pressure test results can be transferred to the InBody.


- * Only compatible with other InBody products
- * Blood pressure values can be printed on the InBody Result Sheet.
- * Serial cable is sold separately.

1. Turn off the BPBIO320S.


* If the BPBIO320S is already turned on while being connected to the InBody, the connection may not be smooth.

2. Connect the serial cable to the RS232C (9-pin serial) connection terminal located on the back of the BPBIO320S. Connect the other end of the serial cable to the serial connection terminal on the InBody device.

* Refer to the InBody manual to find the location of the serial connection terminal.

Example: InBody770

3. Turn on the BPBIO320S.

4. Devices will connect automatically.

Note

- Refer to the InBody product manual to print the blood pressure values on the InBody Result Sheet.

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I. Precautions for Maintenance

Caution

- Never move the equipment while the power is on.
- Do not let foreign objects, like dust and liquid to enter the device. Foreign objects can damage electric components.
- Clean the exterior of the device gently with a lint-free cloth once a week. Do not scratch the LED screen while cleaning the equipment.
- Replace the cuff cover when it becomes dirty or worn out.
- Packing material and other waste should be disposed of according to the laws and regulations.

Check the following before turning the equipment on.

- Make sure the equipment was not damaged from impact.
- Make sure the equipment is not contaminated.
- Make sure the equipment is not wet.
- Make sure the power cord is completely plugged in.
- Make sure the power cord is not damaged.
- Make sure the power cord is not a tripping hazard.

The BPBIO320S calibrates automatically when it is turned on. But it may be necessary for a specialist to inspect the device once every two years to maintain its performance and safety. If an error message appears during the selftest, please contact product support.

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A. Precautionary Steps

Caution


- User should rest at least 5 to 10 minutes before the test.
- Do not roll up the user's sleeve, but take off thick clothes.
- User should sit up straight in the chair.
- Insert the user's arm until the elbow is caught in the elbow groove.
- Adjust the height of the chair so that the inserted upper arm and the heart are horizontal.
- Test in a relaxed state.
- User should not talk or move during the test.
- The BPBIO320S uses the right arm. Test results may be inaccurate if the left arm is used.

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B. Test Instructions


1. Adjust the height of the chair so that the inserted upper arm and the heart are horizontal.


2. Rest for at least 5 minutes prior to taking the test.

* If examined in a bit unstable condition, the test result may be incorrect, and generally it will turn up higher than your usual blood pressure.

3. Insert the arm until the elbow is caught in the elbow groove, as shown.

4. Press  button to start the test.

5. Pressure is automatically applied to the cuff, and the blood pressure is measured.

* Do not talk or move during the test.


* Relax while sitting up straight.

Caution

- If an emergency occurs during the test, press  button or the  button on the front bottom left of the BPBIO320S. The cuff will release quickly.
- If the air is not released after the button is pressed, turn off the BPBIO320S.

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
6. The cuff is automatically released after completing the test. Wait for the cuff to fully loosen and pull out your arm.
7. Check the results on the screen.

- * The voice guide will be heard and the Result Sheet will be printed depending on settings.
- * Do not pull the Result Sheet while it is printing. It will cut automatically.

* Example of test result screen

* Example of Result Sheet


(To print last five blood pressure test results, press and hold the print button while in standby mode.)

PULSE Graph(F12): OFF PULSE Graph(F12): ON PULSE Graph(F12): OFF PULSE Graph(F12): ON
Option(F13): OFF Option(F13): OFF Option(F13): ON Option(F13): ON


Last 5 tests

3 consecutive measurements mode (F12 : ON)

2017/12/06 19:12
SYSTOLIC 110 mmHg
DIASTOLIC 75 mmHg
PULSE 96 bpm

2017/12/06 19:12
SYSTOLIC 110 mmHg
DIASTOLIC 75 mmHg
M.A.P 87 mmHg
P.P 35 mmHg
PULSE 96 bpm
R.P.P 10560

TIME	SYS	DIA	PUL
17/12/06 09:09	111 mmHg	078 mmHg	103 bpm
17/12/06 09:10	116 mmHg	075 mmHg	094 bpm
17/12/06 09:10	118 mmHg	074 mmHg	098 bpm
17/12/06 09:11	117 mmHg	080 mmHg	097 bpm
17/12/06 09:12	110 mmHg	075 mmHg	096 bpm

Average of 3 measurements		
SYSTOLIC	110 mmHg	
DIASTOLIC	75 mmHg	
PULSE	96 bpm	
SYS	DIA	PUL
1. 105	79	96
2. 111	75	96
3. 115	72	96


Note

- Refer to 'Setup' for change Result Sheet settings.

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- * The user can check results by using 'Check Blood Pressure Result' at the bottom right of the Information Panel
- * Consult a physician for accurate diagnosis.

Note


- Rest for 5 minutes before testing again.
- Blood pressure can vary. Consult your physician for accurate diagnosis.

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
C. Test Posture


For correct test results, sit with proper posture while taking the test.

*Be sure to place the elbow in the elbow groove before starting the test.


*Remove thick clothes. Do not roll up sleeves.

III. Transportation and Storage

A. Cautions during Transportation

To transport the BPBIO320S safely, two people should keep the equipment parallel to the ground.





Caution

- When transporting the equipment alone, keep the equipment horizontal while transporting.

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

B. Repacking Instructions

1. Turn off the BPBIO320S.
2. Remove all connected cables and put the packing material into the BPBIO320S. Wrap package in protective vinyl. Insert the protective box into the anti-theft groove on the rear of the BPBIO320S.

3. Place the large air cushion on the bottom of the packing box and place the BPBIO320S on it. Cover with the small air cushion and seal the box with tape.

⚠ Caution

- When repacking the equipment, protective packing materials provided by InBody must be used.

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C. Transportation and Storage Environment

The BPBIO320S should be transported or stored under the following conditions.


Temperatures range	-20 to 70°C (-4 to 158°F)
Relative humidity	10 to 95% RH (No Condensation)
Atmospheric pressure range	50 to 106 kPa

IV.FAQ

If a problem persists, please contact InBody.

A.Error Code

The error messages can be displayed on the screen if they occur during the test.




Please provide the error code to Product Support.

B.Regarding the Equipment

If problems cannot be resolved with the information below, please contact InBody.

Problem	Cause	Corrective measures
The device is not turning on.	This can occur when the power cable is not fully inserted into the power outlet.	<ul style="list-style-type: none"> Plug the power cable into a grounded 3-prong outlet. <div style="display: flex; justify-content: space-around; align-items: flex-start;"> <div style="text-align: center;"> <p><Plug Type B></p> </div> <div style="text-align: center;"> <p><Plug Type F></p> </div> </div> <ul style="list-style-type: none"> When using a power strip, make sure the switch is turned on. <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p><Plug Type B></p> </div> <div style="text-align: center;"> <p><Plug Type F></p> </div> </div> <ul style="list-style-type: none"> Make sure the power cable (AC 250V 10A) provided by InBody is in use. <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p><Plug Type B></p> </div> <div style="text-align: center;"> <p><Plug Type F></p> </div> </div>

Question	Cause	Measures
The device is not turning on.	The fuse is blown.	<p>Turn off the BPBIO320S, and replace the fuse.</p> <p>❶ The fuse socket is located on the back of the device beneath the power plug.</p>  <p>❷ Pull the fuse out by using a flathead screwdriver. Replace the fuse with the provided spare fuse. Additional fuses can be purchased if necessary.</p>  <p>TYPE: Fast-Acting Rated current: 2.5AL Rated voltage: 250V</p> <p>⚠ Warning</p> <ul style="list-style-type: none"> • Be sure to turn off the power when replacing the fuse. • The correct rated fuse must be used. (250V, 2.5AL) <p>⚠ Avertissement</p> <ul style="list-style-type: none"> • Éteignez complètement l'appareil lorsque vous devez remplacer le fusible. • Veuillez à utiliser un fusible de calibre approprié (250V, 2.5 A)
		<p>Turn on the volume knob on the rear of the BPBIO320S.</p> 
Voice guide does not work.	The volume knob on the rear of the blood pressure monitor may be turned off.	
	The sound setting has been set to 'Beep guide' or 'Disable speaker.'	<p>If the sound setting <i>F01</i> is set at <i>b e e f</i>, only the 'beep' sound (alarm sound) will be heard. If it is set at <i>o f f f</i>, no sound will be heard. If you want the voice guide to be heard, set the sound setting to <i>u o i</i>.</p> <p>* Refer to Setup on page 13 to change the sound settings.</p>

Question	Cause	Measures
The Result Sheet is not printing.	Printing may be disabled.	Check whether the settings for automatic printing option is on/off. If the setting is OFF the result sheets will only be printed manually. (Manual printing: The result sheet is printed only when the  button is pressed.) * Refer to page 13 'F. Setup' section and change the automatic printing option.
	The printer may be out of paper.	Since the printer papers are consumables, the replacement cycles of the paper depend on the number of times result sheets have been printed.
	Printer paper may be installed incorrectly.	Open the cover, and check the direction of the printer paper. * Refer to 'Loading and Reloading Printer Paper' on page 9.
	The cover may not be closed completely.	Make sure the central cover is properly closed.

Caution

- The company is not liable for the problems occurred when others repair and inspect the equipment except employees of InBody Co., Ltd.

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Note







- Since error messages and things of detail such as incorrect result sheets can be used as a basis for A/S, please record or keep all of them.

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
C. Regarding the Test

If problems cannot be resolved with the information below, please contact InBody.

Question	Answer
<ul style="list-style-type: none"> • Test results are incorrect. 	<ul style="list-style-type: none"> • The arm and heart are not aligned during the test. Adjust the height of the chair to align the arm and heart. <div style="display: flex; align-items: flex-start; margin-top: 10px;"> <div style="margin-right: 20px;">  </div> <div> <p>Correct testing posture. The height of the arm and the heart is aligned.</p> </div> </div> <div style="display: flex; align-items: flex-start; margin-top: 20px;"> <div style="margin-right: 20px;">  </div> <div> <p>Incorrect testing posture. The arm is higher than the heart. Blood pressure value will be measured lower.</p> </div> </div> <div style="display: flex; align-items: flex-start; margin-top: 20px;"> <div style="margin-right: 20px;">  </div> <div> <p>Incorrect testing posture. The arm is lower than the heart. Blood pressure value will be measured higher.</p> </div> </div>
<ul style="list-style-type: none"> • Inserting the arm too far or not far enough may produce incorrect results. Place elbow in elbow groove for accurate results. 	<div style="display: flex; align-items: center; margin-top: 10px;"> <div style="margin-right: 20px;">  </div> <div>  </div> </div>
<ul style="list-style-type: none"> • Thick clothes or rolled up sleeves may produce incorrect results. 	<div style="display: flex; align-items: center; margin-top: 10px;"> <div style="margin-right: 20px;">  </div> <div style="margin-right: 20px;">  </div> <div>  </div> </div>
<ul style="list-style-type: none"> • The user may have tested without sufficient rest. Take the test after resting at least 5 minutes. • Do not talk or move during the test. 	

* Blood pressure may vary throughout the day. InBody recommends testing at different times of the day.

Question	Answer
<ul style="list-style-type: none"> When should blood pressure be measured? 	<ul style="list-style-type: none"> Blood pressure may vary throughout the day, as shown below.  <p>Therefore, the best time to test is when users can test regularly on a daily basis. It is recommended to test at the same time each day under the same conditions.</p>
<ul style="list-style-type: none"> How frequently can tests be done? 	<ul style="list-style-type: none"> The automatic blood pressure monitor applies pressure to the arm to measure blood pressure. Consecutive testing may put a strain on blood vessel. Rest for about 5 minutes between tests.
<ul style="list-style-type: none"> Is there any case that I should not test? 	<ul style="list-style-type: none"> Automatic blood pressure monitor applies pressure to the arm to measure blood pressure. If you had blood drawn or had an injection given in the either arm, it is recommended to wait at least 1 hour before testing, or measure with another arm.
<ul style="list-style-type: none"> Can results vary depending on the testing environment and temperature? 	<ul style="list-style-type: none"> BPBIO320S can operate at 5 to 40°C (41 to 104°F). However, if the user feels that the temperature is low or high, this feeling makes the blood vessels expand or contract, so it will ultimately affect the blood pressure results. It is recommended to test in a stable condition in an environment where no cold or heat is felt. The user will not be able to perform a natural and stable test if the user is testing in a loud or distracted environment. It is recommended to test in an environment free from noise and disturbances.
<ul style="list-style-type: none"> The blood pressure values are different at home and at the hospital. 	<ul style="list-style-type: none"> The blood pressure may change: <ul style="list-style-type: none"> - Within 1 hour after meals - After having alcohol or caffeine - After smoking - After taking a bath - After exercise - After urination or defecation - When testing in an unfamiliar location
<ul style="list-style-type: none"> The blood pressure values are different at home and at the hospital. 	<ul style="list-style-type: none"> Blood pressure can change depending on the environment. Test in the same environment when managing blood pressure. Psychologically, people tend to feel stable at home than at hospital. Hospital-tested blood pressure tends to be higher than home-tested blood pressure.


V. Others

- * The BPBIO320S was manufactured by the quality control procedures of InBody Co., Ltd. InBody Co., Ltd. complies with ISO9001 and ISO13485 which are international quality management systems.
- * This equipment complies with IEC60601-1, the international safety standard for electronic medical equipment. This equipment also complies with IEC60601-1-2, the international standard for electromagnetic compatibility.

A.Exterior and Functions

The names and functions of each part of the BPBIO320S are as follows:

- * Please check each component of the system for damage prior to installation.

- | | | |
|----|--------------|--|
| 1 | | : A display screen showing results, current time, etc. |
| 2 | | : A button to print Result Sheet |
| 3 | | : A button to start or stop a test |
| 4 | | : Emergency stop button |
| 5 | | : Printer paper cover |
| 6 | | : Motion sensor to automatically switch from sleep to standby mode |
| 7 | Elbow Sensor | : A sensor to guide the correct examination posture |
| 8 | | : UP/DOWN buttons to adjust settings |
| 9 | | : Volume control knob |
| 10 | | : RS232C (9-pin serial) connection terminal |
| 11 | | : A port used to connect a mercury monitor in order to adjust the pressure value |
| 12 | | : A terminal used to connect a multi-display device |
| 13 | | : Input for external power source (AC 100-240V, 50/60Hz), Fuse (F2.5AL, 250V) |
| 14 | | : A switch used to turn the equipment on/off |
| 15 | | : A groove that can be fixed to a table for theft prevention |
| 16 | | : Speaker |

B. Safety Information

Warning

- If this equipment has been altered, it must be inspected to ensure that it can be used safely.
- Results can only be interpreted by experienced healthcare professionals. It cannot be used for diagnosis, medication, or other treatments performed at the consumer's discretion without a doctor's prescription.

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Caution

- Excessively high or low temperatures, humidity and pressure may affect the operation of the equipment and may cause malfunction. Use blood pressure monitor within the range specified in the product specifications.
- Do not touch the external device connections such as RS232C (9-pin serial) connection terminal on the back of the equipment.
- Packing material and other wastes should be disposed of according to the relevant laws and regulations.

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
Safety Symbols

	Danger high voltage
	Warning / Caution
	Note
	BF-type equipment
	Power on
	Power off

Other Symbols

	Manufacturer		Serial number
	Authorized representative in the EUROPEAN COMMUNITY		Alternating current
	European conformity		Operating instructions
	Complies with AAMI ES60601-1 CSA-C22.2 No. 60601-1 (MET mark only applies for USA and Canada.)		

C. Product Classification

Oscillometric multi-functional electronic blood pressure monitor

- Type of protection against electric shock: Class 1 device
- Level of protection against electric shock: BF-type
- Level of protection against flooding: General Equipment (No special protection against external water infiltration)
- Class against electromagnetic radiation intensity: Class A

D. Specifications

Test Instructions	Oscillometric
Test range	Pressure: 0-300 mmHg Pulse: 30-240 bpm
Degree of Precision	Pressure: ± 2 mmHg Pulse: $\pm 1.5\%$
Test result	Systolic blood pressure, diastolic blood pressure, Mean arterial pressure, Pulse pressure, PULSE rate, RPP (rate pressure product)
Test Time	Approx. 30 seconds on average (20-50 seconds depending on the pulse and blood pressure value)
Pressurization Time	Approx. 10 seconds
Minimum scale unit	1 mmHg
Screen Type	7-Segment LED
Result output	Result value, PULSE graph, Last 5 measurements
Storage function	Can store results for up to 1 million tests
Energy saving	Automatically enters the energy saving mode after 2 minutes of idle time.
Dual safety mechanism	When the Start/Stop button is pressed, the cuff rapidly releases air. The cuff automatically releases when the air pressure exceeds 300 mmHg. When the Emergency Stop button is pressed, the cuff quickly deflates.
Voice guide	Guides user through blood pressure measurement and results.
Printer	High-speed thermal printer with a built-in automatic cutter (width: 2.25 in / 5.7 cm)
Rated voltage and power consumption	AC 100-240V, 50-60Hz, 30VA
Dimensions	489 (W) \times 409 (D) \times 284 (H) mm 19.3 (W) \times 16.1 (L) \times 11.2 (H) in
Weight	9kg (19.84 lb)
Operating environment	Temperature 5 \sim 40°C (41 \sim 104°F), Humidity 15 \sim 90% RH, Atmospheric pressure 70 \sim 106 kPa
Storage environment	Temperature -20 \sim 70°C (-4 \sim 158°F), Humidity 10 \sim 95% RH, Atmospheric pressure 50 \sim 106 kPa (No Condensation)

* Specifications are subject to change without notice.

* This product is a medical device. Please read the guide and precautions before using this blood pressure monitor.

E. EMC Declaration

Test Name	IEC 60601-1-2 4th Edition	Ref. Standard	Ports to Test	AC Mains Voltage	Test Level Required	Notes
Mains terminal disturbance voltage	4th	CISPR 11:2015/ AMD2:2019	AC Mains	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 230V, 50Hz	Group 1, Class B	
Radiated disturbance	4th	CISPR 11:2015/ AMD2:2019	Enclosure	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 230V, 50Hz Battery	Group 1, Class B	
Harmonic Current Emissions	4th	IEC 61000-3-2:2018	AC Mains	230V~50Hz	Class A	
Voltage Fluctuations and Flicker	4th	IEC 61000-3-3:2013/ AMD1:2017	AC Mains	230V~50Hz	Pst = 1 Plt = 0.65 Dmax = 4% DC = 3.3%	
Electrostatic Discharges (ESD)	4th	IEC 61000-4-2:2008	Enclosure	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 230V, 50Hz Battery	±8 kv contact ±2 kv, ±4 kv ±8 kv, ±15 kv Air	1 pulse/ 1sec contact 8 kV air 15kV
Radiated RF electromagnetic Fields	4th	IEC 61000-4-2:2008	Enclosure	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 230V, 50Hz Battery	10 V/m 80 MHz-2.7 GHz 80% AM at 1kHz RF Wireless Comm. (Refer to test report clause 1.15)	Dwell Time is 3 sec
Electric Fast Transients and bursts	4th	IEC 61000-4-4:2012	AC Mains	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 230V, 50Hz	±2 kv, AC, 100 hk PRR	5 pulses at 0°, 90°, 180°, 270°
Surges	4th	IEC 61000-4-5:2014/ AMD1:2017	AC Mains	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 230V, 50Hz	±0.5 kv, ±1 kv L to N (CM) ±0.5 kv, ±1 kv ±2 kv L & N to PE (CM)	5 pulses at 0°, 90°, 180°, 270°
Conducted Disturbances, induced by RF fields	4th	IEC 61000-4-6:2013	AC Mains & all I/O	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 230V, 50Hz Battery	3 Vrms 150 Hk to 80 Mk 80% AM at 1Hk, 2Hz 6 Vrms in ISM Bands, amateur radio bands between 0.15 Mk and 80Mk	Dwell time is 3 sec

Power-frequency Magnetic Field	4th	IEC 61000-4-8:2009	Enclosure	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 230V, 50Hz Battery	30 A/m 50 Hz & 60 Hz	
Voltage Dips, Interruptions, and variations	4th	IEC61000-4-11:2004/ AMD1:2017	AC Mains	100V, 50Hz 100V, 60Hz 120V, 60Hz 220V, 60Hz 240V, 50Hz 240V, 60Hz	0% UT for 0.5 cycle	AT 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°
					0% UT for 1 cycle	AT 0°, 180°
					50Hz : 70% UT for 25 cycles 60Hz : 70% UT for 30 cycles	AT 0°, 180°
					50Hz : 70% UT for 250 cycles 60Hz : 70% UT for 300 cycles	AT 0°, 180°

Electromagnetic immunity

The [Model name] is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the [Model name]. Otherwise, degradation of the performance of this equipment could result.

Immunity test	Band a)	Service a)	Modulation b)	IEC60601 test level	Compliance level
Proximity fields from RF wireless Communications IEC61000-4-3	380 - 390 MHz	TETRA 400	Pulse modulation b) 18Hz	27 V/m	27 V/m
	430 – 470 MHz	GMRS 460 FRS 460	FM c) ±5 kHz deviation 1 kHz sine	28 V/m	28 V/m
	704 – 787 MHz	LTE Band13, 17	Pulse modulation b) 217 Hz	9 V/m	9 V/m
	800 – 960 MHz	GSM800:900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation b) 18 Hz	28 V/m	28 V/m
	1700 – 1990 MHz	GSM 1800 CDMA 1900 GSM 1900 DECT LTE Band 1,2,4,25 UMTS	Pulse modulation b) 217 Hz	28 V/m	28 V/m
	2400 – 2570 MHz	Bluetooth WLAN 802.11b/g/n RFID 2450 LTE Band 7	Pulse modulation b) 217 Hz	28 V/m	28 V/m
	5100 – 5800 MHz	WLAN 802.11a/n	Pulse modulation b) 217 Hz	9 V/m	9 V/m

NOTE : If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1m. The 1m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50% duty cycle square wave signal.

c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.